

HOW TO TRIM YOUR EARMOULD TUBING



•Step 1: Your ear-mould will arrive with long tubing, like this

•Step 2: Put the ear-mould in your ear, ensuring you insert it correctly.



•Step 3: Put the hearing aid behind the ear and mark the point where the hearing aid meets the tubing (see arrow)



•Step 4: Trim the tubing about quarter of an inch/3mm up from the marked point and then connect the tubing to the hearing aid elbow (hook)

Department contact information:

Phone: 02921 843179

Email: Audiology.Helpline.CAV@wales.nhs.uk

Address: Audiology Clinic, Clinic 9, Outpatients Department, University Hospital of Wales, Heath Park Way, Cardiff, CF14 4XW

Walk in repair clinic is currently closed please post repairs to the above address