



GIG
CYMRU
NHS
WALES

Bwrdd Iechyd Prifysgol
Caerdydd a'r Fro
Cardiff and Vale
University Health Board

UNIVERSITY HOSPITAL LLANDOUGH

ALL WALES

ADULT CYSTIC FIBROSIS CENTRE

How to keep your nebuliser equipment clean

- 1) Fill the base with 80 ml of cold water



- 2) Load your **washed** nebuliser equipment into machine



- 3) Switch on to start cycle



Wait 5 minutes after the steriliser has finished before emptying content, the steriliser will remain hot for sometime