

**APOMORPHINE ADMINISTERING PROCEDURE VIA SUB-CUTANEOUS
INFUSION DEVICE TO IN-PATIENTS WITH PARKINSON'S DISEASE**

Reference No:	UHB 184	Version No:	1	Previous Trust / LHB Ref No:	T/276
----------------------	---------	--------------------	---	-------------------------------------	-------

Documents to read alongside this Procedure	N/A
---	-----

Classification of document:	Corporate
Area for Circulation:	In-patient clinical areas
Author/Reviewee:	Parkinson's Disease Specialist Nurse / Medical Equipment Training Officer
Executive Lead:	Medical Director
Group Consulted Via/ Committee:	Medicines Management Group
Approved by:	Medicines Management Group
Date of Approval:	25 February 2013
Date of Review:	25 February 2015
Date Published:	26 April 2013

Disclaimer

When using this document please ensure that the version you are using is the most up to date either by checking on the UHB database for any new versions. If the review date has passed please contact the author.

OUT OF DATE POLICY DOCUMENTS MUST NOT BE RELIED ON

Version Number	Date of Review Approved	Date Published	Summary of Amendments
1	25/02/2013	26/04/2013	Document updated. Supersedes previous Trust procedure reference no 276

PROCEDURE ADMINISTERING APOMORPHINE VIA SUB-CUTANEOUS INFUSION DEVICE TO IN-PATIENTS WITH PARKINSON'S DISEASE

- Patients who are receiving Apomorphine therapy will normally use a specific ambulatory device, the APO-GO pump, which is set up by their carers; however staff in the acute hospital setting are not trained to use this device so an alternative, the McKinley T34, is to be used.
- Please see overleaf for the full procedure/checklist involved in setting up a McKinley T34 syringe driver to administer Apomorphine in an acute hospital setting.
- Do not alter this procedure unless advised to do so by the Parkinson's disease Nurse or Doctor.
- Apomorphine stains green – these stains are resistant to most cleaning methods, be careful during set up and disposal of syringes.
- To reduce skin irritation it is important to alternate the site of injection daily. Three sites are used – the anterior abdominal wall and the anterior aspect of each thigh.

If you have any queries or experience any problems regarding the guidance in this document, please contact:-

The Parkinson's Disease Nurse Specialist can be contacted on Rookwood Ext. 33838

Clinical Engineering Help Desk (UHW Ext. 45678) for information on the T34 Syringe Driver.

Pharmacy (UHW Ext. 42988; Llandough Ext. 25261) for information on Apomorphine.

CHECKLIST FOR SETTING UP A MCKINLEY T34 SYRINGE DRIVER FOR APOMORPHINE

Place tick in box

- | | | |
|-----|--|--------------------------|
| 1. | Apomorphine is available from pharmacy in pre-filled 10ml syringes (5mg/ml). | <input type="checkbox"/> |
| 2. | Obtain a McKinley T34 and disposables pack from the Medical Equipment Loan Service. | <input type="checkbox"/> |
| 3. | <p>Calculate the correct volume of Apomorphine for the patient's prescription:
(The usual dosage is between 4-7mg/hr)</p> <p style="margin-left: 200px;">For example.
4.5mg/hr over 15 hours.
4.5 x 15 = 67.5mg = 13.5ml
Therefore draw up 14.5ml to compensate for priming the line.</p> | <input type="checkbox"/> |
| 4. | <p>Calculate the rate required (in this case after priming):</p> $rate = \frac{volume\ (ml)\ remaining\ in\ syringe}{time\ in\ hours}$ <p>Draw up the required dose of Apomorphine, from the pre-filled syringes:</p> | <input type="checkbox"/> |
| 5. | <p>-Take the Apo-go connector that is provided with the pre-filled syringe from Pharmacy and detach one of the caps. Screw open end onto the luer lock syringe.</p> <p>-Detach the second connector cap. Remove the grey rubber cap from the Pre-Filled Syringe and attach to the connector.</p> <p>-Holding the syringes vertically, push down on the plunger of the pre-filled syringe to transfer the required volume of Apomorphine.</p> | <input type="checkbox"/> |
| 6. | Prepare a label and attach to the syringe. | <input type="checkbox"/> |
| 7. | Examine T34 to ensure it is serviceable. | <input type="checkbox"/> |
| 8. | Insert battery – switch on T34. | <input type="checkbox"/> |
| 9. | <p>Wait for self check to complete:</p> <p style="margin-left: 20px;">Display lights up</p> <p style="margin-left: 20px;">Alarm sounds</p> <p style="margin-left: 20px;">Drive mechanism moves to last position</p> <p style="margin-left: 20px;">Configured data is displayed</p> | <input type="checkbox"/> |
| 10. | Connect the butterfly set to syringe and prime the line | <input type="checkbox"/> |
| 11. | Position the drive mechanism on the T34 to fit syringe length using the ◀◀ and the ▶▶ keys. | <input type="checkbox"/> |
| 12. | Fit the syringe and confirm the type and size. | <input type="checkbox"/> |
| 13. | Insert the butterfly device subcutaneously and secure with a transparent dressing. | <input type="checkbox"/> |
| 14. | Confirm the volume measured by the T34 (this should be the volume to be infused) | <input type="checkbox"/> |
| 15. | Enter the calculated rate in ml/hr using the ▲ and ▼ keys. | <input type="checkbox"/> |
| 16. | Confirm the rate entered. | <input type="checkbox"/> |
| 17. | <p>The display will show the current infusion data:</p> <p style="margin-left: 20px;">Volume xx.xml</p> <p style="margin-left: 20px;">Duration xx:xx</p> <p style="margin-left: 20px;">Rate x.xml/hr</p> | <input type="checkbox"/> |
| 18. | If correct confirm the settings. | <input type="checkbox"/> |
| 19. | Start the infusion & lock the keypad [L]. | <input type="checkbox"/> |
| 20. | Place the syringe driver <u>correctly</u> into its box and lock. | <input type="checkbox"/> |
| 21. | Place the syringe driver out of direct sunlight at approximately the same height as the infusion site. Document date, time and volume left in syringe. | <input type="checkbox"/> |
| 22. | <p>Check syringe driver after 10 minutes to:</p> <ul style="list-style-type: none"> - ensure that the volume is decreasing appropriately - the site is not inflamed - that the battery capacity is greater than 30% | <input type="checkbox"/> |